

The Lipo Cannulae Catalogue



 gpsurgicalsurgicalservices.com

Welcome

Leading-Edge Expertise

Our representatives are trained in advanced techniques and latest generations of equipment. Known for their trusted expertise, leadership and dedication to bringing the highest quality service and supplies to our surgical customers. Our valued customers are leaders in their areas of expertise, they collaborate with us for extraordinary care and improved patient outcomes.

We Believe in Building Great Relationships

Attention to detail is imperative in the operating room.

We understand. That's why our highest priority is getting to know you and your operating room team. Your personalities, your working styles, your patients. We believe this sets the foundation for a long-term relationship which creates value-added benefits that not only include access to quality supplies, training and resources from our highly-trained team, but also creates operational efficiencies in your practice.

Our aim is to supply the NHS and Private healthcare providers with a Quality, Reliable and Cost Effective service.

We are always striving to achieve the highest possible standards of quality surgical instruments, this is why we work very closely with the end users to achieve our goals.

Tanvir Akhtar

Tanvir Akhtar
Sales Director

GP Surgical Services Limited
Unit 6 Honeycrest Industrial Park
Staplehurst
Kent
TN12 0RXX



Table of Content

Liposuction Cannula Handles	1-2
Cannula Hole Pattern	3-6
Infiltration Cannulas	7
Reinjection Micro Cannula	8-9

Liposuction Cannula Handles



BLACK HANDLE

Product Code: GPSS-2201

Product Description:

Black Handle

Handle made of aluminum with type 2 black anodized, 6.65" length 1" diameter



BLUE HANDLE

Product Code: GPSS-2202

Product Description:

BLUE Handle

Handle made of aluminum with type 2 blue anodized, 6.65" length 1" diameter



GOLD HANDLE

Product Code: GPSS-2203

Product Description:

Gold Handle

Handle made of aluminum with type 2 gold anodized, 6.65" length 1" diameter



HANDLE w/SWIVEL BARB

Product Code: GPSS-2204

Product Code: GPSS-2205

Product Description:

Handle made of aluminum with type 2 blue anodized, 6.65" length, 1" diameter
(for Large Handle, 3/4" for Small Handle)



HANDLE

Product Code: GPSS-2206

Product Code: GPSS-2207

Product Description:

Handle made of aluminum with type 2 blue anodized, 6.65" length 1" diameter
(for Large Handle, 3/4" for Small Handle)

Liposuction Cannula Handles (cont.)



STAINLESS STEEL TOOMEY 60cc HUB

Product Code: GPSS-2208

Product Description:
Stainless Steel and Stainless Steel hub



ALUMINUM TOOMEY 60cc HUB

Product Code: GPSS-2209

Product Description:
Hub made of Aluminum with type 2 blue anodized, 1.55" length 1.31" diameter



QUICK SCREW HUB

Product Code: GPSS-2210

Product Description:
Hub made of Stainless Steel, 2.1" length 7/8" diameter



SMALL QUICK DISCONNECT HUB

Product Code: GPSS-2211

Product Description:
Hub made of Stainless Steel, 1.93" length 5/8" diameter



QUICK DISCONNECT HUB

Product Code: GPSS-2212

Product Description:
Hub made of Stainless Steel, 1.6" length 7/8" diameter



QUICK LOCK HUB

Product Code: GPSS-2213

Product Description:
Hub made of Stainless Steel, 1.65" length 1/2" diameter



LUER LOCK HUB

Product Code: GPSS-2214

Product Description:
Luer Lock Hub
Hub made of Stainless Steel, 1.16" length 3/8" diameter

Cannula Hole Pattern



MERCEDES

Product Code: GPSS-2215
Available Diameter & Length
Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.
Length: 15cm, 23cm, 26cm, 30cm, 32cm.



TRI PORT

Product Code: GPSS-2216
Available Diameter & Length
Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.
Length: 15cm, 23cm, 26cm, 30cm, 32cm.



LAS VEGAS

Product Code: GPSS-2217
Available Diameter & Length
Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.
Length: 15cm, 23cm, 26cm, 30cm, 32cm.



STANDARD 3 PORT

Product Code: GPSS-2218
Available Diameter & Length
Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.
Length: 15cm, 23cm, 26cm, 30cm, 32cm.



STANDARD 2 PORT

Product Code: GPSS-2219
Available Diameter & Length
Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.
Length: 15cm, 23cm, 26cm, 30cm, 32cm.

Cannula Hole Pattern (cont.)



STANDARD SINGLE PORT

Product Code: GPSS-2220

Available Diameter & Length

Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.

Length: 15cm, 23cm, 26cm, 30cm, 32cm.



3 PORT 60 DEGREE INDEX

Product Code: GPSS-2221

Available Diameter & Length

Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.

Length: 15cm, 23cm, 26cm, 30cm, 32cm.



3 PORT 60 DEGREE INDEX

Product Code: GPSS-2222

Available Diameter & Length

Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.

Length: 15cm, 23cm, 26cm, 30cm, 32cm.



BASKET

Product Code: GPSS-2223

Available Diameter & Length

Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.

Length: 15cm, 23cm, 26cm, 30cm, 32cm.



SPATULA

Product Code: GPSS-2224

Available Diameter & Length

Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.

Length: 15cm, 23cm, 26cm, 30cm, 32cm.

Cannula Hole Pattern (cont.)



SPATULA 2 PORT

Product Code: GPSS-2225

Available Diameter & Length

Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.

Length: 15cm, 23cm, 26cm, 30cm, 32cm.



SPATULA 3 PORT

Product Code: GPSS-2226

Available Diameter & Length

Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.

Length: 15cm, 23cm, 26cm, 30cm, 32cm.



SPATULA 3 PORT

Product Code: GPSS-2227

Available Diameter & Length

Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.

Length: 15cm, 23cm, 26cm, 30cm, 32cm.



COBRA

Product Code: GPSS-2228

Available Diameter & Length

Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.

Length: 15cm, 23cm, 26cm, 30cm, 32cm.



FUENTES

Product Code: GPSS-2229

Available Diameter & Length

Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.

Length: 15cm, 23cm, 26cm, 30cm, 32cm.

Cannula Hole Pattern (cont.)



SAYLAN

Product Code: GPSS-2230

Available Diameter & Length

Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.

Length: 15cm, 23cm, 26cm, 30cm, 32cm.



CANDY CANE

Product Code: GPSS-2231

Available Diameter & Length

Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.

Length: 15cm, 23cm, 26cm, 30cm, 32cm.



MODERN SPIRAL

Product Code: GPSS-2232

Available Diameter & Length

Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.

Length: 15cm, 23cm, 26cm, 30cm, 32cm.



TOLEDO SINGLE PORT

Product Code: GPSS-2233

Available Diameter & Length

Diameter: 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm, 5mm.

Length: 15cm, 23cm, 26cm, 30cm, 32cm.

Infiltration Cannulas



GARDEN SPRAY ANESTHESIA NEEDLE

Product Code: GPSS-2234

Available Diameter & Length

Diameter: 18ga, 17ga, 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm.

Length: 15cm, 23cm, 26cm, 30cm, 32cm.



MULTI HOLE SPRAY NEEDLE

Product Code: GPSS-2235

Available Diameter & Length

Diameter: 18ga, 17ga, 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm.

Length: 15cm, 23cm, 26cm, 30cm, 32cm.



BODY ANESTHESIA NEEDLE

Product Code: GPSS-2236

Available Diameter & Length

Diameter: 18ga, 17ga, 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm.

Length: 15cm, 23cm, 26cm, 30cm, 32cm.



Infusion 6 hole Needle

Product Code: GPSS-2237

Available Diameter & Length

Diameter: 18ga, 17ga, 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm.

Length: 15cm, 23cm, 26cm, 30cm, 32cm.



Infusion 8 hole Needle

Product Code: GPSS-2238

Available Diameter & Length

Diameter: 18ga, 17ga, 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm.

Length: 15cm, 23cm, 26cm, 30cm, 32cm.

Reinjection Micro Cannula



MICRO INJECTION NEEDLE, STYLE 3

Product Code: GPSS-2239

Available Diameter & Length

Diameter: 18ga, 17ga, 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm.

Length: 15cm, 23cm, 26cm.



MICRO INJECTION NEEDLE, STYLE 2

Product Code: GPSS-2240

Available Diameter & Length

Diameter: 18ga, 17ga, 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm.

Length: 15cm, 23cm, 26cm.



MICRO INJECTION NEEDLE, STYLE 1

Product Code: GPSS-2241

Available Diameter & Length

Diameter: 18ga, 17ga, 16ga, 15ga, 14ga, 13ga, 12ga, 3mm, 4mm.

Length: 15cm, 23cm, 26cm, 30cm, 32cm.

Reinjection Micro Cannulas (cont.)



MICRO INJECTION NEEDLE "V" DISSECTOR

Product Code: GPSS-2242

Available Diameter & Length
Diameter: 18ga, 17ga, 16ga.
Length: 5cm, 7cm, 8cm, 9cm, 12cm.



BEVEL 45

Product Code: GPSS-2243

Available Diameter & Length
Diameter: 16ga, 15ga, 14ga.
Length: 7cm, 9cm, 12cm, 15cm, 23cm.



BEVEL 30

Product Code: GPSS-2244

Available Diameter & Length
Diameter: 16ga, 15ga, 14ga.
Length: 7cm, 9cm, 12cm, 15cm, 23cm.

terms and conditions of sale

"The Company" refers to "gpsurgicalservices.com" "the Customer" means the individual, practice, firm, company, or other party with whom the Company contracts.

This agreement is applicable to all goods and a service supplied in the United Kingdom by "Company", whether ordered from this catalogue or otherwise and supersede or cancels all previous printed Terms and Conditions of Sale.

1. CUSTOMER CONDITIONS.

These conditions apply to any contract of sale between the Company and the Customer to the exclusion of any different terms and conditions issued by the Customer, except insofar as the Company and the Customer may have otherwise agreed.

gpsurgicalservices.com does not supply goods to the general public via its website. However, enquiries are invited from the general public to info@gpsurgicalservices.com

2. PRICES.

While every effort is made to maintain prices, the Company reserves the right to alter prices without notice at any time before delivery, to reflect increases imposed by manufacturers and when currency exchange rates adversely affect the cost of imported goods.

3. V.A.T.

Prices are exclusive of V.A.T. at 20% (where applicable). The Company will vary prices to reflect any changes in V.A.T. rate of V.A.T. due in respect of the sale.

4. CARRIAGE.

Carriage and packing will be charged on orders under £50 at £7.95 per box of less than 25kgs (excluding V.A.T.).

5. PAYMENT.

Full payment must be made within 30 days of statement date. The Company reserves the right to remove discounts for late payments.

6. OVERDUE ACCOUNTS.

Orders received from a customer whose account is overdue may not be processed until the account is settled in full.

7. DELIVERY.

Any date or period for delivery or any date of delivery stated in any contract is intended by the Customer and accepted by the Buyer as being an estimate only, not giving rise to contractual obligations.

8. SHORTAGES, DAMAGES AND NON DELIVERY.

Should goods not arrive within 14 days of the despatch dates shown on the invoice, the Company must be advised immediately so that a claim may be made on the carrier. In the case of damage or shortages the Company must be notified within 3 working days of receipt of goods. In case of a claim all cartons and packing material should be retained for examination either by the carrier or insurance company.

9. CANCELLATION OF GOODS.

The Company shall endeavour to meet Customer & individual requirement but the Company shall be under no obligation to accept cancellation or amendment of any order or any part of an order. Where such cancellation or amendment is accepted it is on the understanding that a cancellation or amendment charge may be levied.

10. PROOF OF DELIVERY.

Where a customer requests such proof and this is provided, the Company reserves the right to charge for any costs incurred in obtaining this information.

11. GOODS RETURN PROCEDURE.

Subject to prior arrangement goods may be accepted back for credit or replacement providing the goods were not obtained specially and that they are in sound and saleable condition. The name and address of the sender must be enclosed with each package and accompanied by an advice note showing full details. A handling charge could be incurred for late returns.

12. SPECIFICATIONS.

Descriptions and illustrations are accurate at time of publication. Due to continuing development, the Company reserves the right to supply later versions of the goods described. Every effort has been made to ensure the accuracy of the information contained in this catalogue. However, the Company cannot be held responsible for any mistakes or for any damages resulting from the use of materials described herein.

13. CONDITIONS AND WARRANTIES.

The following provisions apply in contracts for the sale of goods as defined by Section 62 of the Sale of Goods Act 1893 as amended hereinafter referred to as "The Act." (a) In the case of a consumer sale as defined by section 55 (7) of the Act, the goods are supplied with the benefit where the same are applicable of the conditions and warranties implied by Sections 12-15 inclusive of the Act. (b) In any other case goods and materials are supplied with the benefit (where the same are applicable) of the conditions and warranties implied by Section 12 of the Act. (c) Save as aforesaid and save as herein expressly set out the goods are not supplied with, or subject to, any conditions warranty or other term whether express or implied unless specifically stated by the Company in writing.

14. RETENTION OF TITLE.

Goods that are supplied by the Company to the Buyer shall remain the property of the Company until full payment for the goods supplied to the Buyer and all other payments whatsoever due from the Buyer to the Company have been made. The separate identity of these goods should be maintained by the Buyer. The Buyer may, however, sell such goods in the normal course of business as agent and bailee for the Company and the proceeds of the sale shall be held separately to the account of the Company up to the amount of the full purchase price of the goods sold. The Customer's right to sell the goods shall terminate immediately upon whichever is the earlier of, written notice from the Company or the appointment of any Receiver, Manager or Administrative Receiver to the whole or any part of the undertaking, business or assets of the Customer or if any petition is presented applying to the Court for the winding up of, or for an Administration Order to be made in relation to, the Customer or a resolution is passed for its winding up or dissolution, or the Customer becomes bankrupt or enters into a composition or voluntary arrangement with its creditors.



WEB SITE ADDRESS:- GPSURGICALSERVICES.COM

INFO@GPSURGICALSERVICES.COM

TELEPHONE: 01622 828854

GPSURGICALSERVICES LTD

UNIT 6 HONEYCREST INDUSTRIAL PARK, LODGE ROAD,
STAPLEHURST, TN12 0RX

